

AIDAL® PLUS

STERILANT
ARTG NO: 135545

AIDAL® PLUS is suitable for use in Hospitals, Dental and Medical Surgeries and Veterinary Hospitals, for the Sterilisation and High Level Disinfection of Critical and Semi Critical Medical Devices.

AIDAL® PLUS is a polyglycol dialdehyde complex instrument disinfectant intended for general use in a surgical practice where sterilisation or high level disinfection is required. It is a stable, odour suppressed, mildly acidic glutaraldehyde sterilant and instrument disinfectant (high level), which exhibits potent biocidal action.

AIDAL® PLUS shows several major improvements over alkaline activated Glutaraldehyde. Improved biocidal performance is achieved with reduced glutaraldehyde odour and reduced propensity for irritation and yellowing to human skin compared with alkaline glutaraldehyde preparations.

AIDAL® PLUS has been tested in a wide array of tests and various conditions, including testing at the end of the shelf life, at the end of the re-use period, diluted to the Minimum Recommended Concentration (MRC) of active material (1.5% Glutaraldehyde), and then soiled and tested in hard water against hardest to kill organisms.

AIDAL® PLUS will not damage endoscopes or other complex medical equipment when used in accordance with directions and the instrument manufacturers instructions.

AIDAL® PLUS has an MRC of 1.5% w/w glutaraldehyde. All tests for sterilant or disinfectant efficacy were conducted at the MRC or lower in accordance with strict performance standards.

BENEFITS:

- AIDAL® PLUS has demonstrated performance as a bactericidal, fungicidal, virucidal, mycobactericidal and sporicidal high level disinfectant with complete efficacy under these conditions as a Class IIb Medical Device.
- AIDAL® PLUS TEST STRIPS to monitor the concentration of glutaraldehyde.

DIRECTIONS FOR USE:

HIGH LEVEL DISINFECTION

Soak instruments in undiluted solution.

35° C for 5 minutes.

25° C for 20 minutes.

Discard the solution immediately if any of the following are noted – the solution changes colour or becomes cloudy, the solution is subject to copious dilution or if the solution is compromised in any way. Use ADIAL PLUS Test Strips to verify the in use solution is above the MRC.

STERILISATION

Soak instruments in undiluted solution.

35° C for 7 hours 40 minutes.

EFFICACY DATA AND SOAK TIMES

AIDAL® PLUS is effective against HIV and HBV.

Due to the many varied test conditions and regulatory requirements in different countries please consult the Front Label of AIDAL® PLUS for the relevant soak time and temperature in the country of use.

DISINFECTION PROCEDURES

Thorough precleaning with an acceptable product is an essential element of any effective cleaning and disinfecting protocol. An appropriate procedure will involve at least 5 steps – preclean, rinse, disinfect, sterile rinse and dry, then store. For further technical advice please contact your local Whiteley representative.

It is of critical importance that only recommended detergents be used in conjunction with AIDAL® PLUS. International authorities recommend that only a neutral detergent, a neutral enzyme cleaner or a biofilm remover be used as part of thorough mechanical cleaning.

MATRIX® has been shown in internationally peer reviewed published clinical data to provide excellent soil and biofilm removal from materials used in endoscopes and other medical devices.

MEDIZYME is a neutral enzyme cleaner which is also suitable for cleaning medical devices. Separate technical literature is available on each product and should be studied prior to use.

It is also of critical importance that all instruments be thoroughly rinsed prior to immersion in AIDAL® PLUS. Care should be taken to ensure that instruments are as dry as possible prior to immersion. AIDAL® PLUS should be carefully drained from the instruments following an appropriate immersion time. The final rinse should be with sterile water using at least 5 times the volume required to cover the instrument.

This rinsing procedure is best done using enough water to immerse the instrument each time.

PRECAUTION WITH METAL, PLASTIC & RUBBER INSTRUMENTS

AIDAL® PLUS solutions should not be used with carbon steel instruments since they may produce a mild corrosion (as does distilled water) after a few hours. 24 hour continuous exposure to solutions will not affect most plastics, brass, aluminium, chrome, nickel or stainless steel. Chrome and nickel plated instruments will not be tarnished after 24 hours, with the exception of nickel plating over a beryllium base. Poorly plated, worn or corroded instruments with cracks could obviously be tarnished after processing.

AIDAL® PLUS is not suitable for use with products made of natural rubber or plasticised PVC.

AIDAL® PLUS contains a corrosion inhibitor; however, it is recommended that disinfected instruments do not remain in the AIDAL® PLUS solution for longer than required for disinfection.

The immersion time for optical instruments should not exceed the soaking time that the equipment manufacturer recommends.

Regarding the general matter of instruments tarnishing, the manufacturer does not warrant anything beyond the above statements since it is well known that the quality of the metals, alloys, plastic parts or components can greatly vary among instrument brands. When in doubt, the customer should contact the instrument manufacturer.

SOLUTION LIFE

Manual soaking or use in automated endoscopic reprocessors is recommended.

AIDAL® PLUS is chemically stable for 15 months, when stored in the original container.

AIDAL® PLUS is intended for re-use up to 28 days, or until indicated by AIDAL® PLUS test strips, or until the solution is visibly soiled or becomes cloudy or discoloured, whichever comes first. AIDAL® PLUS should not be reused beyond 28 days even if test strips indicate a pass result.

AIDAL® PLUS may be reused for up to 28 days for high level disinfection procedures (after opening or putting in an AER). In re-use situation AIDAL® PLUS test strips should be used to test the solution.

AIDAL® PLUS TEST STRIPS

Whiteley recommend the use of AIDAL® PLUS TEST STRIPS to monitor the concentration of glutaraldehyde. During re-use it is logical that rinse water introduced into AIDAL® PLUS will cause the concentration of glutaraldehyde to fall, therefore it is advisable to test the concentration at a minimum of each day or as recommended by your Gastroenterological Association. If the concentration of AIDAL® PLUS falls below the MRC of 1.5% the AIDAL® PLUS TEST

STRIPS will give a negative reading. If there is any doubt over the result another test strip should be used to verify the result. If there is any remaining doubt over re-use of AIDAL® PLUS the solution should be discarded immediately. AIDAL® PLUS TEST STRIPS are available from your AIDAL® PLUS supplier.

STORAGE AND DISPOSAL

This product must be kept in locked storage to make it inaccessible for children or persons unfamiliar with its proper use. AIDAL® PLUS should be stored at or below 25°C.

AIDAL® PLUS contains only biodegradable materials. Glutaraldehyde has been shown to be without effect at concentrations below 10 ppm. Where possible, disposal should be with copious quantities of water. Where septic systems are involved the active glutaraldehyde may be neutralised with dilute solutions of sodium bisulphite before disposal. Tests reveal that glutaraldehyde is approximately 90% biodegradable within a standard five day biodegradable test at levels below microbiological effectiveness.

The container may be disposed of by thoroughly rinsing with water and then discarding through the normal waste system. It is recommended that the container be punctured prior to disposal to prevent re-use.

The bottle, lid and label are all made from fully recyclable materials where recycling is available.

Where required, residual glutaraldehyde can be neutralised by adding neutraliser. Special instructions apply.

ACTIVE CONSTITUENT

21 g/L glutaraldehyde in the form of buffered neutral 1, 5-pentanedial poly glycol complex.

WARNING

May cause sensitisation by skin contact or inhalation.

AIDAL® PLUS is hazardous according to the criteria of Safework Australia. Glutaraldehyde can act as an irritant and sensitiser to the skin and irritates the airways of the respiratory tract. Appropriate ventilation, work practices and personal protective equipment are necessary to prevent skin contact and inhalation.

Occupational asthma and/or rhinitis have been indicated in a number of workers exposed to glutaraldehyde and other aldehyde solutions.

Carcinogenicity: Glutaraldehyde is not a carcinogen.

Reproductive Toxicity / Teratogenicity: There is no difference of incidence of miscarriage between those exposed to glutaraldehyde and those not exposed to the chemical.

Mutagenicity: The results of more recent assays have generally shown that glutaraldehyde is mutagenic in vitro. All in vivo tests to date have been negative. Consequently, glutaraldehyde does not meet the criteria for classification as a mutagen.

SAFE WORKING RECOMMENDATIONS

Glutaraldehyde is known to cause allergic reactions and skin sensitisation through skin contact or inhalation.

Occupational asthma and/or rhinitis have been indicated in a number of workers exposed to glutaraldehyde.

Use only in covered airtight containers to avoid evaporation and contamination. Avoid contact with eyes, skin and food.

Gloves should be worn at all times [Double Gloving or Nitrile gloves are internationally recommended]. Protective clothing eg. Long sleeved shirts and plastic aprons and covered impervious shoes should be worn when using glutaraldehyde.

Avoid breathing vapour. The area in which the product is used should be well ventilated at all times. Glutaraldehyde is a respiratory irritant.

The Safety Data Sheet on AIDAL® PLUS provides additional safety data, and should be read before using AIDAL® PLUS. Consult the Material Safety Data Sheet for first aid information.

ABOUT LMS:

Our comprehensive learning platform offers a broad range of information including the latest information on Infection Prevention and Control, recommended products and cleaning processes for individual markets.



Scan to register for Whiteley's
Online Learning Centre



SAFETY DATA: Information is listed in the Safety Data Sheet.

PACK SIZE: 2 x 5 Litre

ENVIRONMENTAL STATEMENT:

Whiteley manufactures products that comply with Australian and International environmental policy. Whiteley aims to ensure that every product is 100% biodegradable, whilst achieving optimal performance and benefit to customers. All Whiteley products are supplied in 100% recyclable containers. Whiteley and its staff are fully committed to ensuring that our environment and those international environs in which our products are used are enhanced by the use of technology, and that these environments are preserved for the enjoyment of future generations.

© 2024 This document remains the property of Whiteley Corporation Pty Ltd.

Alterations are not permitted without prior written authorisation from Whiteley Corporation Pty Ltd.